

WonderWink® Women's Consultation Lab Coat WW4072



- 65/35 poly/cotton poplin
- Locker loop
- Left chest pocket with pen slot
- Top-loading front pockets
- Interior tablet pocket
- Back yoke
- 3-button placket opening
- Center back vent
- WonderWink heat transfer logo at left pocket
- Size medium: 29.5-inch back length

++

CARE INSTRUCTIONS


Machine Wash Cold With Like Colors. Gentle Cycle. Use Only Non-Chlorine Bleach When Needed. Tumble Dry Low. Remove Promptly From Dryer. Cool Iron When Needed.



fron



HOW TO MEASURE



BUST

Measure under the arm and around the fullest part of the bust with arms down, keeping tape horizontal.

SIZE CHART

| | 2XS | XS | S | М | L | XL | 2XL | 3XL | 4XL |
|------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| Size | 0 | 2 | 4-6 | 8-10 | 12-14 | 16-18 | 20-22 | 24-26 | 28-30 |
| Bust | 31-32 | 33-34 | 35-36 | 37-39 | 40-43 | 44-47 | 48-51 | 52-55 | 56-59 |

COLOR INFORMATION

White PMS WHITE